

THE MORTARA ELI10

12 lead rest electrocardiograph



Power Supply:	100-240 V 50/60 Hz, 50 V Internally rechargeable battery
Handheld unit:	is charged via the included docking station
Instrument Type	12-lead Electrocardiograph
Input Channels	Simultaneous acquisition of 12 leads
Standard Leads Acquired	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
Waveform Display	Backlit 1/4 VGA LCD, 3-channel lead groups or complete 8-channel presentation
Input Impedance ANS/AAMI	Meets or exceeds the requirements of EC11 Input Dynamic Range Electrode Offset Tolerance
ECG Storage	Up to 60 ECGs
Dimensions	portable unit 17,0 x 12,8 x 3,5 cm
Weight	4.9 kg including battery
Dimensions portable ECG	17,0 x 12,8 x 3,5 cm

Prices

Price will be based on the agreed minimum rental period and on specifications of required equipment.

Our standard minimum rental terms are 6, 12, 18, 24, 36, 48 and 60 months.

For a suitable tailor made quotation, please contact us at:

- E: sales@medicapitalrent.com
- T: +31 71 581 4 581

Should you require different technical specifications, please do not hesitate to contact us.

Easy to use

The Mortara Eli 10 is an easy-to-use reliable 12-lead diagnostic electrocardiograph capable of viewing, acquiring, printing, and storing ECG data. Ideal for use in medical practices and small laboratories.

Backlit LCD Display

High-resolution display provides real time preview of 12-lead ECG and post-acquisition review of acquired ECG.

ECG Printing

ECG printing is enabled through a cardiology management system using the network or a stand-alone USB external printer.

Connectivity Options

Multiple connectivity options standard LAN and USB port in docking station. Available extra options include internal modem, GSM/GPRS mobile or WLAN.

Internal Storage

Internal storage for up to 60 digital ECG records. Via docking station or optional external USB memory.

ECG quality control

This cardiograph provides extensive ECG quality-control features including AC noise filtering. The display allows preview of any three leads at a time to verify ECG quality before recording. The patient cables are well insulated to insure safety.

MediCapital Rent BV | Tieleman en Drosweg 17 - 19 | 2314 XW | Leiden
The Netherlands | T. +31 (0)71 581 45 81 | F. +31 (0)71 581 4 582

MediCapital Rent INC | 101 TDK Boulevard, Suite F Peachtree City | GA 30269
United States | T. +1 (678) 969 22 07 | F. +1 (770) 248 05 00

info@medicapitalrent.com | www.medicapitalrent.com



MediCapital Rent
Your reliable partner in clinical research