

THE MORTARA ELI250

12 lead rest electrocardiograph



Power Supply: 100-240 V 50/60 Hz, 50 V Internally rechargeable battery

Instrument Type: 12-lead Electrocardiograph

Input Channels: Simultaneous acquisition of 12 leads

Standard Leads: Acquired I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
Waveform Display Backlit 1/4 VGA LCD, 3-channel lead groups or complete 8-channel presentation

Input Impedance: Meet or exceed the requirements of **input**

Dynamic Range: ANSI/AAMI EC11

Electrode Offset Tolerance

Common Mode Rejection

ECG Storage: Up to 60 ECGs

Weight: 5.1 kg including battery

Dimensions: 39.4 x 43.2 x 10.2 cm

Prices

Standard pricing is Flat Rate (no minimum rental period required).

Discounted pricing available on the basis of agreed minimum rental period upfront.

Available minimum rental terms are 6, 12, 18, 24, 36, 48 and 60 months.

For a suitable tailor made quotation, please contact us at:

- E: sales@medicapitalrent.com
- T: +31 71 581 4 581

Should you require different technical specifications, please do not hesitate to contact us.

Easy-to-use

The Mortara Eli 250 is an easy-to-use reliable 12-lead diagnostic electrocardiograph capable of viewing, acquiring, printing, and storing ECG data. Ideal for use in medical practices and small laboratories.

Backlit LCD Display

High-resolution display provides real time preview of 12-lead ECG and post-acquisition review of acquired ECG.

Full Keypad

Alphanumeric elastomer keypad also includes dedicated "onetouch" buttons for ECG acquisition, rhythm printing and ECG transmission.

Full-size, Z-fold Paper

Full-size paper (8.5" x 11" or A4 supported) output includes multiple Standard or Cabrera print options including 3+1, 3+3, 6-channel and 12-channel outputs.

Internal Storage

Internal storage for up to 60 digital ECG records. Optional expanded storage for up to 150 digital ECG records.

ECG quality control

This cardiograph provides extensive ECG quality-control features including AC noise filtering. The display allows preview of any three leads at a time to verify ECG quality before recording. The patient cables are well isolated to insure safety.

MediCapital Rent BV | Tieleman en Drosweg 17 - 19 | 2314 XW | Leiden
The Netherlands | T. +31 (0)71 581 45 81 | F. +31 (0)71 581 4 582

MediCapital Rent INC | 101 TDK Boulevard, Suite F Peachtree City | GA 30269
United States | T. +1 (678) 969 22 07 | F. +1 (770) 248 05 00

info@medicapitalrent.com | www.medicapitalrent.com



MediCapital Rent

Your reliable partner in clinical research